

---

---

फिनाइल जे-अम्ल, तकनीकी — विशिष्टि  
( दूसरा पुनरीक्षण )

**Phenyl J-Acid,  
Technical — Specification**  
( *Second Revision* )

ICS 71.080.40

© BIS 2021




भारतीय मानक ब्यूरो  
BUREAU OF INDIAN STANDARDS  
मानक भवन, 9 बहादुरशाह ज़फर मार्ग, नई दिल्ली – 110002  
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG  
NEW DELHI-110002  
[www.bis.gov.in](http://www.bis.gov.in) [www.standardsbis.in](http://www.standardsbis.in)

## FOREWORD

This Indian Standard (Second Revision) was adopted by the Bureau of Indian Standards after the draft finalized by Dye Intermediates Sectional Committee had been approved by the Petroleum Coal and Related Products Division Council.

Phenyl J-acid, chemically described as 2-phenglamino-5-naphthol-7-sulphonic acid ( $C_{16}H_{13}O_4NS$ ) is an important intermediate used in the manufacture of dyes. It is also known as 2-phenylamino-5-hydroxy naphthalene-7-sulphonic acid. It is represented by the following structural formula:



N-PHENYL-J-ACID  
(MOLECULAR MASS 315 )  
CAS Registry Number [119-40-4]

This standard was first published in 1975 and subsequently revised in 1994. In the first revision, the requirements of matter insoluble in dilute sodium carbonate and J-acid content were modified.

Considering to the development in latest analytical techniques in the last 15 years, committee decided to revise this standard to incorporate purity by high performance liquid chromatography and impurities of aniline and J acid content in N-phenyl J acid by HPLC method.

The composition of the Committee, responsible for the formulation of this standard is given at Annex D.

For the purpose of deciding whether a particular requirement of this standard is complied with the final value, observed or calculated, expressing the result of a test or analysis, shall be rounded off-in accordance with IS 2 : 1960 'Rules for rounding off numerical values (*revised*)'. The number of significant places retained in the rounded off value should be the same as that of the specified value in this standard.

# Indian Standard

## PHENYL J-ACID, TECHNICAL — SPECIFICATION

( Second Revision )

### 1 SCOPE

This standard prescribes the requirements and the methods of sampling and test for phenyl J-acid, technical.

### 2 REFERENCES

The following Indian Standards contain provisions which through reference in the text, constitute provisions of this standard. At the time of publication the editions indicated were valid. All standards are subject to revision, and parties to agreement, based on the standard are encouraged to investigate the possibility of applying the most recent editions of the standards indicated below:

<i>IS No.</i>	<i>Title</i>
797 : 1982	Common salt for chemical industries ( <i>third revision</i> )
1070 : 1992	Reagent grade water ( <i>third revision</i> )
2552 : 1989	Steel drums (galvanized and ungalvanized) — Specification ( <i>third revision</i> )
5299 : 2001	Methods of sampling and tests for dye intermediates ( <i>first revision</i> )

### 3 REQUIREMENTS

#### 3.1 Description

The material shall be in the form of a paste, or, if dry, in the form of light grey to greenish grey powder.

**3.2** The material shall also comply with the requirements given in Table 1.

### 4 PACKING AND MARKING

#### 4.1 Packing

The material shall be packed in steel drums (see IS 2552) lined with suitable polyethylene film, or as agreed to between the purchaser and the supplier.

**4.2** Each container shall be securely closed and shall bear legibly and indelibly the following information:

- Name of the material;
- Indication of the source of manufacture;

- Net mass of material;
- Month and year of the manufacture;
- Batch or lot number; and
- The minimum cautionary notice worded.

“IT IS A MILD SENSITIZER LOCAL CONTACT  
MAY CAUSE DERMATITIS”

**Table 1 Requirements for Phenyl  
J-Acid, Technical**

( *Clauses 3.2, 5.3.1, 5.3.2 and 6.1* )

Sl No.	Characteristic	Requirement	Method of Test, Ref to
(1)	(2)	(3)	(4)
i)	Purity by Coupling Value, <i>Min</i>	85	A
ii)	Purity by HPLC (on dry basis)	85	B
iii)	Aniline content, percent by mass (HPLC method), <i>Max</i>	0.2	C
iv)	J-acid content, percent by mass (HPLC method), <i>Max</i>	0.5	C
v)	Matter insoluble in dilute sodium carbonate solution, percent by mass, <i>Max</i>	0.3	11.2 of IS 5299

**4.2.1** The containers may also be marked with the Standard Mark.

#### 4.2.2 BIS Certification Marking

The product(s) conforming to the requirements of this standard may be certified as per the conformity assessment schemes under the provisions of the *Bureau of Indian Standards Act, 2016* and the Rules and Regulations framed thereunder, and the products may be marked with the Standard Mark.

### 5 SAMPLING

**5.1** Representative samples of the material shall be drawn as prescribed in 4 of IS 5299.

#### 5.2 Numbers of Tests

**5.2.1** Tests for assay shall be conducted on each of the individual samples separately.

**5.2.2** Tests for the determination of the remaining characteristics, namely, aniline content, J-acid content and matter insoluble in dilute sodium carbonate solution shall be conducted on the composite sample.

### **5.3 Criteria for Conformity**

#### **5.3.1 For Individual Samples**

The lot shall be declared as conforming to the requirement of assay if each of the individual test results satisfies the relevant requirement given in Table 1.

#### **5.3.2 For Composite Sample**

The lot shall be declared as conforming to the requirements of aniline content, J-acid content and

matter insoluble in dilute sodium carbonate solution if the test results satisfy the relevant requirements given in Table 1.

## **6 TEST METHODS**

**6.1** Tests shall be conducted according to methods prescribed in Annexes and IS Specifications as given in column 4 of Table 1.

### **6.2 Quality of Reagents**

Unless specified otherwise, 'pure chemicals' and distilled water (*see* IS 1070) shall be employed in tests.

NOTE — 'Pure chemicals' shall mean chemicals that do not contain impurities which affect the results of analysis.

## ANNEX A

[ Table 1, SI No. (i) ]

## DETERMINATION OF PURITY BY COUPLING VALUE METHOD

## A-0 OUTLINE OF METHOD

The material is dissolved in dilute sodium carbonate solution. A known volume of the solution is coupled with standard 4-nitroaniline diazo in sodium carbonate and dilute sodium acetate medium and from the consumption.

## A-1 PREPARED SAMPLE

Dry the material at  $105 \pm 1$  °C to constant mass. Grind and mix well. Transfer the material to a wide-mouthed bottle and stopper it. Do not expose the sample to an atmosphere containing acidic or alkaline fumes. Use this prepared sample for tests.

## A-2 REAGENTS

## A-2.1 Alkaline H-Acid Indicator Solution

## A-2.2 Sodium Acetate Solution — 20 percent (m/v).

## A-2.3 Sodium Carbonate Solution — 10 percent (m/v).

## A-2.4 4-nitroaniline Diazo Solution — 0.1 N.

## A-2.5 Sodium Nitrite 1.0 N Solution

## A-2.6 Common Salt — see IS 797.

## A-2.7 Preparation of the PNA-Diazo Solution 0.1 mol/l (0.1 N)

6.9066 g of *p*-nitroaniline are weighed into a 400 ml glass beaker and dissolved in 100 ml of purify water and 50 ml of hydrochloric acid 30-32 percent. The solution is warmed and stirred until all compounds are completely dissolved. After, the solution is cooled down with ice to 3-5 °C. At this point, 50.00 ml of sodium nitrite, concentration ( $\text{NaNO}_2$ ) = 1.0 mol/l are added directly into the solution in order to avoid a loss of  $\text{NO}_2$  ( $\text{NO}_x$  vapours.). The solution is poured into a 500 ml measuring flask and taken to volume with purify water. Keep the solution in ice bath.

A-2.8 Standardization of 0.1 N PNA Diazo Solution:-  
Weight 1.0 g (accuracy: 0.0001 g)  $\beta$ -naphthol and transfer in 500 ml beaker, add 10-15 ml of 1N

NaOH to dissolved. Stir well if not dissolved. Add 200 ml distilled water . After addition of 30 ml sodium acetate solution, and add ice to cool it ,the solution is titrated gradually to the end point with *p*-nitro aniline-diazo (PNA-diazo) . Test with alkaline H-acid for excess diazo and with diazo for the coupling component. A faint test line with alkaline H-acid reproducible for a period of 5 min denotes the end point.

## A-2.9 Calculation

Normality of PNA-diazo solution: -

$$\frac{99 \times \text{weight of } \beta\text{-Naphthol [STD] (in gm)}}{V \times M}$$

where

V = volume in ml, of the PNA-diazo solution used;  
and

M = mass in g, of  $\beta$ -naphthol taken for the test.

## A-3 PROCEDURE

Weight accurately about 1 g of the prepared sample (see A-1) and transfer to a 500 ml beaker with the help of water. Add 20 ml of sodium carbonate solution and 100 ml of sodium acetate solution. Stir mechanically and cool with ice to 0 to 5 °C. Titrate with standard 4 nitroaniline diazo solution, from a cold water-jacketted burette. Add 15 g of common salt near the end point. Test with alkaline H-acid for excess diazo and with diazo for the coupling component. A faint test line with alkaline H-acid reproducible for a period of 5 min denotes the end point.

## A-4 CALCULATION

$$\text{Assay, percent by mass, on dry basis} = \frac{V \times N \times 31.5}{M}$$

where

V = Volume in ml, of the standard 4 nitroaniline diazo solution used;

N = Normality of 4 nitroaniline diazo solution; and

M = Mass in g, of the prepared sample taken for the test.

## ANNEX B

[ Table 1, Sl No. (ii) ]

## PURITY BY HPLC

**B-0 OUTLINE OF METHOD**

High-performance liquid chromatography or high-pressure liquid chromatography (HPLC) is a chromatographic method that is used to separate a mixture of compounds in analytical chemistry and biochemistry so as to identify, quantify or purify the individual components of the mixture.

**B-1 OBJECTIVE**

To determine purity of phenyl J-acid by high performance liquid chromatography.

**B-2 APPARATUS**

Binary gradient liquid chromatography system with UV detector capable of being operated under conditions suitable for resolving the individual constituents into distinct peak may be used.

**B-3 COLUMN**

Kinetex C18 100A 250 × 4.6 mm, 5 µm or equivalent.

**B-4 REAGENTS****B-4.1 Acetonitrile, HPLC Grade****B-4.2 Water, HPLC Grade****B-4.3 Di-sodium Hydrogen Orthophosphate, HPLC Grade****B-4.4 Ammonium Di-hydrogen Orthophosphate, HPLC Grade****B-4.5 Phenyl J-acid, Known Purity****B-5 STANDARD PREPARATION**

Weigh accurately 0.0500 gm standard phenyl J-acid in 100 ml volumetric flask dissolve it in water: acetonitrile (20:80) and make up to the mark with water: acetonitrile (20:80).

**B-6 SAMPLE PREPARATION**

Weigh accurately 0.0500 gm sample in 100 ml volumetric flask dissolve it in water: acetonitrile (20:80) and make up to the mark with water: acetonitrile (20:80).

**B-7 BUFFER PREPARATION**

Take 0.5750 gm ammonium di-hydrogen orthophosphate and 0.7000 gm di-sodium hydrogen orthophosphate in 1 liter volumetric flask. Add 200 ml HPLC grade water and complete dissolve it. Make total volume with HPLC grade water.

**B-8 FLOW RATE: 1.00 ml/min****B-9 MOBILE PHASE:**

Acetonitrile	Buffer
25	75

**B-10 COLUMN OVEN TEMPERATURE: 40 °C****B-11 INJECTION VOLUME: 5 µl****B-12 RUN TIME: 15 min****B-13 WAVE LENGTH: 230 nm**

<b>B-14 PEAK TIME:</b>	J-acid	- 2.34 min
	Aniline	- 5.81 min
	Phenyl J-acid	- 8.35 min

**B-15 CALCULATION**

Calculate the peak area of individual constituent pertaining to phenyl J-acid on the chromatogram of the material. The concentration of the constituent may be obtained on the basis peak area on chromatogram obtained with known amount of pure phenyl J-acid.

$$\text{Percent of phenyl J-acid} = \frac{A2 \times V1 \times W1 \times B2}{A1 \times V2 \times W2 \times B1} \times 100$$

where

A1	=	Area of standard phenyl J-acid;
V1	=	Injection volume of standard phenyl J-acid;
W1	=	Weight of standard phenyl J-acid;
B1	=	Total volume of standard phenyl J-acid;
A2	=	Area of phenyl J-acid peak in sample;
V2	=	Injection volume of sample;
W2	=	Weight of sample;
B2	=	Total volume of sample;




FIG. 1

## ANNEX C

[ Table 1, Sl No. ( iii) and (iv) ]

### DETERMINATION OF J-ACID AND ANILINE CONTENT BY HPLC

#### C-1 OBJECTIVE

To determine J-acid and aniline content in N-phenyl J-acid by high performance liquid chromatography.

#### C-2 APPARATUS

Binary gradient liquid chromatography system capable of being operated under conditions suitable for resolving the individual constituents into distinct peak may be used.

#### C-3 COLUMN

C18 100A 250 × 4.6 mm, 5 µm or equivalent.

#### C-4 REAGENT

##### C-4.1 Acetonitrile, HPLC Grade

##### C-4.2 Water, HPLC Grade

##### C-4.3 Di-sodium Hydrogen Orthophosphate, HPLC Grade

##### C-4.4 Ammonium Di-hydrogen Orthophosphate, HPLC Grade

##### C-4.5 Phenyl J-acid, Known Purity

#### C-5 STANDARD PREPARATION

**C-5.1** Weigh accurately 0.0050 gm standard J-acid in 100 ml volumetric flask dissolve it in water: acetonitrile (20:80) and make up to the mark with water: acetonitrile (20:80). ----- A

**C-5.2** Take 1 ml of solution A dilute 10 ml with water : acetonitrile (20:80).

**C-5.3** Weigh accurately 0.0050 gm standards aniline in 100 ml volumetric flask dissolve it in acetonitrile and make up to the mark with acetonitrile. ----- B

**C-5.4** Take 1 ml of solution B dilute 10 ml with acetonitrile.

#### C-6 SAMPLE PREPARATION

Weigh accurately 0.0500 gm sample in 100 ml volumetric flask dissolve it in water: acetonitrile (20:80) and make up to the mark with water: acetonitrile (20:80).

#### C-7 BUFFER PREPARATION

Take 0.5750 gm ammonium di-hydrogen orthophosphate and 0.7000 gm di-sodium hydrogen orthophosphate in 1 litre volumetric flask. Add 200 ml HPLC grade water and complete dissolve it. Make total volume with HPLC grade water.

**C-8 FLOW RATE:** 1.00 ml/min

#### C-9 MOBILE PHASE:

Acetonitrile	Buffer
25	75

**C-10 COLUMN OVEN TEMPERATURE:** 40 °C

**C-11 INJECTION VOLUME:** 5 µl




FIG. 2

**C-12 RUN TIME:** 15 min

**C-13 WAVE LENGTH:** 230 nm

**C-14 PEAK TIME:**

J Acid - 34 min  
 Aniline - 5.81 min  
 Phenyl J Acid - 8.35 min

#### C-15 CALCULATION

Calculate the peak area of individual constituent pertaining to J-acid/aniline on the chromatogram of the material. The concentration of the constituent may be obtained on the basis peak area on chromatogram obtained with standard J-acid/aniline.

Percent of J-acid/aniline =

$$\frac{A2 \times V1 \times W1 \times B2}{A1 \times V2 \times W2 \times B1} \times 100$$

where

A1 = Area of standard J-acid/aniline;  
 V1 = Injection volume of standard J-acid/aniline;  
 W1 = Weight of standard J-acid/aniline;  
 B1 = Total volume of standard J-acid/aniline;  
 A2 = Area of Jacid/aniline peak in sample;  
 V2 = Injection volume of sample;  
 W2 = Weight of sample; and  
 B2 = Total volume of sample.



**ANNEX D***( Foreword )***COMMITTEE COMPOSITION**

Dye Intermediates Sectional Committee PCD 26

<i>Organization</i>	<i>Representative(s)</i>
Personal Capacity	PROF P. M. BHATE ( <b>Chairman</b> )
Aarti Industries Limited, Mumbai	SHRI KIRIT H. DESAI
Alkyl Amines Chemicals Limited, Mumbai	SHRI KIRAT PATEL SHRI S. V. NIKUMBHE ( <i>Alternate</i> )
Ankleshwar Research and Analytical Infrastructure Limited (ARAIL), Ankleshwar, Gujrat	SHRI MANSUKH H. VEKARIA
Archroma, Thane	DR RAJESH RAMAMURTHY SHRI HEMANT MHADESHWAR ( <i>Alternate</i> )
Atul Limited, Atul, Valsad, Gujarat	DR J. G. DESAI DR M. U. RAHMAN ( <i>Alternate</i> )
BASF India Limited, Mumbai	SHRI UDAY KULKARNI
Central Revenue Control Laboratory (CRCL), New Delhi	DR T. A. SREENIVASA RAO SHRI PRAFUL DALAL ( <i>Alternate</i> )
Colourtex Industries Limited, Mumbai	DR PANKAJ DESAI SHRI ISMAIL HABIBULLA ( <i>Alternate</i> )
Deepak Nitrite Limited, Vadodara	SHRI SAILASH RAVAL DR MANGALYA KAR ( <i>Alternate I</i> ) DR J. K. ASTIK ( <i>Alternate II</i> )
Ecological and Toxicological Association of Dyes (ETAD), Vadodra	DR PARITI SIVA RAMA KUMAR
Gujarat Dyestuffs Manufacturers Association (GDMA), Ahmedabad	SHRI R. S. PATEL SHRI BANSIBHAI PATAL ( <i>Alternate</i> )
Gujarat Pollution Control Board (GPCB), Gandhinagar, Ahmedabad	SHRI Y. A. TAI
Heubach Colour Private Limited, Mumbai	SHRI J. I. SEVAK SHRI VINOD MADHUKAR ( <i>Alternate</i> )
Indian Chemical Council, Mumbai	SHRI P. S. SINGH
Jay Chemical Industries, Ahmedabad	SHRI VILPESH YADAV SHRIMATI MAITRI VYAS ( <i>Alternate</i> )
Ministry of Chemicals & Fertilizers, Department of Chemicals & Petrochemicals, New Delhi	SHRI JASBIR SINGH
Ministry of Defence (Research and Development), New Delhi	DR VIJAY TAK
National Test House (er), Kolkata	SHRI P. K. CHAKRABORTY SHRI Y. C. NIJAWAN ( <i>Alternate</i> )
Sudarshan Chemical Industries Limited, Pune	DR JAYCHANDRAN M. NAIR

<i>Organization</i>	<i>Representative(s)</i>
The Bombay Textile Research Association	DR ASHOK N. DESAI
The Dyestuffs Manufacturers Association of India (DMAI), Mumbai	SHRI V. R. KANETKAR
BIS Director General	SHRI N. K. BANSAL, SCIENTIST 'F' AND HEAD (PCD) [REPRESENTING DIRECTOR GENERAL ( <i>Ex-officio</i> )]

*Member Secretary*  
SHRI CHANDRAKESH SINGH  
SCIENTIST 'D' (PCD), BIS



## Bureau of Indian Standards

BIS is a statutory institution established under the *Bureau of Indian Standards Act, 2016* to promote harmonious development of the activities of standardization, marking and quality certification of goods and attending to connected matters in the country.

### Copyright

BIS has the copyright of all its publications. No part of these publications may be reproduced in any form without the prior permission in writing of BIS. This does not preclude the free use, in the course of implementing the standard, of necessary details, such as symbols and sizes, type or grade designations. Enquiries relating to copyright be addressed to the Director (Publications), BIS.

### Review of Indian Standards

Amendments are issued to standards as the need arises on the basis of comments. Standards are also reviewed periodically; a standard along with amendments is reaffirmed when such review indicates that no changes are needed; if the review indicates that changes are needed, it is taken up for revision. Users of Indian Standards should ascertain that they are in possession of the latest amendments or edition by referring to the latest issue of 'BIS Catalogue' and 'Standards: Monthly Additions'.

This Indian Standard has been developed from Doc No.: PCD 26 (14042).

### Amendments Issued Since Publication

Amend No.	Date of Issue	Text Affected

## BUREAU OF INDIAN STANDARDS

### Headquarters:

Manak Bhavan, 9 Bahadur Shah Zafar Marg, New Delhi 110002

Telephones: 2323 0131, 2323 3375, 2323 9402

Website: [www.bis.gov.in](http://www.bis.gov.in)

### Regional Offices:

Telephones

Central	: Manak Bhavan, 9 Bahadur Shah Zafar Marg NEW DELHI 110002	{ 2323 7617 2323 3841
Eastern	: 1/14 C.I.T. Scheme VII M, V.I.P. Road, Kankurgachi KOLKATA 700054	{ 2337 8499, 2337 8561 2337 8626, 2337 9120
Northern	: Plot No. 4-A, Sector 27-B, Madhya Marg CHANDIGARH 160019	{ 265 0206 265 0290
Southern	: C.I.T. Campus, IV Cross Road, CHENNAI 600113	{ 2254 1216, 2254 1442 2254 2519, 2254 2315
Western	: Manakalaya, E9 MIDC, Marol, Andheri (East) MUMBAI 400093	{ 2832 9295, 2832 7858 2832 7891, 2832 7892

**Branches :** AHMEDABAD. BENGALURU. BHOPAL. BHUBANESHWAR. COIMBATORE.  
DEHRADUN. DURGAPUR. FARIDABAD. GHAZIABAD. GUWAHATI.  
HYDERABAD. JAIPUR. JAMMU. JAMSHEDPUR. KOCHI. LUCKNOW.  
NAGPUR. PARWANOO. PATNA. PUNE. RAIPUR. RAJKOT. VISAKHAPATNAM.

Published by BIS, New Delhi